

# **myGaze® EyeMouse Play User Guide**

**Version 1.5**

January 2015

**visual  
interaction**

## Contents

Part 1	Welcome.....	2
	1.1 Contacting Visual Interaction GmbH.....	3
	1.2 Product Liability.....	4
	1.3 Magnet Precautions.....	5
	1.4 Declaration of Conformity.....	8
	1.5 Product Maintenance.....	9
	1.6 Document Conventions.....	10
	1.7 Glossary.....	11
Part 2	myGaze Assistive System.....	15
	2.1 Minimum System Requirements.....	16
	2.2 Components.....	18
	2.3 Optimal User Conditions.....	18
Part 3	Getting Started.....	21
Part 4	Installing EyeMouse Play.....	25
	4.1 Downloading the Installer.....	25
	4.2 Installing the Software.....	26
Part 5	Mounting the myGaze Assistive Eye Tracker.....	29
	5.1 Reference Point.....	30
	5.2 Mounting on a Desktop Monitor.....	31
	5.3 Mounting on a Laptop Display.....	35
Part 6	EyeMouse Play Overview.....	40

---

6.1	myGaze Menu.....	40
6.2	Access Button.....	42
6.3	Status Colors.....	43
6.4	Gaze Positioning Icon.....	44
6.5	Positioning Window.....	44
6.6	Overriding EyeMouse Play with PC Mouse.....	46
6.7	Updating the Software.....	46
Part 7	Using EyeMouse Play.....	50
7.1	Unhiding the Gaze Positioning Icon.....	50
7.2	Starting EyeMouse Play .....	51
7.3	Positioning Users.....	54
7.4	Opening the myGaze Menu.....	56
7.5	Selecting an Item from myGaze Menu.....	58
7.6	Activating Cursor Mode.....	58
7.7	Pausing EyeMouse Play.....	60
7.8	Activating Dwell-Gaze Mouse Click Mode.....	62
7.9	Activating Double-Click Mode.....	64
7.10	Activating Click and Hold Mode.....	65
7.11	Shutting Down EyeMouse Play.....	67
Part 8	Configuring EyeMouse Play.....	69
8.1	Viewing Configuration Window.....	70
8.2	Managing Profiles.....	71
8.2.1	Creating a Profile.....	71
8.2.2	Selecting a Profile.....	73

---

	8.2.3	Deleting a Profile.....	74
	8.3	Performing a Calibration.....	75
	8.3.1	Calibration Methods.....	76
	8.3.2	Cascading Points.....	78
	8.3.3	Setting Calibration Properties.....	78
	8.3.4	Running a Calibration.....	83
	8.3.5	Improving a Calibration.....	87
	8.4	Setting Mouse Properties.....	88
	8.5	Setting Access Properties.....	91
	8.6	Setting Keyboard Shortcuts.....	93
	8.7	System Settings.....	95
Part 9		EyeMouse Play and Third-Party Software.....	98
	9.1	Grid 2 AAC Overlay.....	98
Part 10		Appendix.....	104
	10.1	Troubleshooting.....	104
	10.2	License Agreement and Warranty.....	105
	10.3	Third Party Licenses.....	108
Part 11		About Visual Interaction.....	110
		Index.....	111

**Welcome**

**Chapter**

**1**

# 1. Welcome

The *myGaze® Assistive System* consists of *myGaze® EyeMouse Play* and *myGaze® Assistive Eye Tracker*, both developed by *Visual Interaction GmbH*. With a few easy steps, this system can be setup and running to assist users with varying abilities to interact with Windows-based applications using only their eye gaze.

This User Guide describes how to install and configure *myGaze® EyeMouse Play* and to mount the *myGaze Assistive Eye Tracker*.



For more information about myGaze®, refer to the *myGaze® User Manual*.

## Document Information

Document Version: 1.5

Release Date: January 2015

Software Version: 1.4.9

For more information as well as access to additional support information and downloads, visit the myGaze website at [www.mygaze.com](http://www.mygaze.com)

Copyright © 2015 Visual Interaction (VI). All other product names are copyright of their respective owners.

See [License Agreement and Warranty](#) for rights and responsibilities of the use of this product.



Please read this User Guide carefully to ensure best results.

## 1.1 Contacting Visual Interaction GmbH

### **Visual Interaction GmbH (VI)**

Behlertstraße 3a/Haus B2  
D-14467 Potsdam  
Germany

Phone: +49 (331) 235 21 52

Fax: +49 (331) 235 22 22

e-mail: [info@mygaze.com](mailto:info@mygaze.com)

Web: [www.mygaze.com](http://www.mygaze.com)

Technical Support: [support@mygaze.com](mailto:support@mygaze.com)

## **1.2 Product Liability**

*Visual Interaction GmbH (VI)* does not assume liability for resultant damages to property or personal injury if the product has been misused in any way or damaged by improper use or failure to observe these operating instructions. In addition, any unauthorized modifications or repairs of the device will render the warranty null and void!



## 1.3 Magnet Precautions

The myGaze Assistive Eye Tracker contains Neodymium magnets (Rare Earth magnets). This allows the device to be quickly connected to the Mounting Bracket which is attached to a monitor. However, this type of magnet is extremely strong and must, therefore, be handled with extreme care.



### Handling Warnings

- Do not confuse Neodymium magnets with standard "fridge magnets". Neodymium magnets can cause injury if not used properly.
- Do not place your fingers between the two magnets while connecting the myGaze Assistive Eye Tracker to the Mounting Bracket. Two attracting magnets have enormous strength and can severely pinch your fingers if placed between the magnets while connecting the myGaze Assistive Eye Tracker.
- Do not let the myGaze Assistive Eye Tracker freely connect to the Mounting Bracket. Although Neodymium magnets have high strength, they are also very brittle and prone to cracking and chipping. If connected too quickly, or if the myGaze Assistive Eye Tracker is dropped, the magnets may shatter and scatter shards of the magnet, possibly even towards the user's eyes.
- Neodymium magnets have strong magnetic fields and likely to cause damage to magnetic media devices. Therefore, keep the myGaze Assistive Eye Tracker away from magnetic media such as hard drives, memory sticks, credit cards, magnetic I.D. cards, or other magnetic media.

- KEEP THE DEVICE IN THE STORAGE CASE WHEN NOT IN USE.
- Do not place metal items near the magnets on the myGaze Assistive Eye Tracker. Metal items such as keys, knives, or tools may cause the magnet to shatter.
- Do not leave the myGaze Assistive Eye Tracker near an open flame or a heat source. Not only will the device be destroyed, but the Neodymium magnets will ignite, burn and create toxic fumes.

## Health Warnings



Neodymium magnets should NEVER be used near a person who uses medical aids such as a pacemaker. The magnet can cause the medical aid to malfunction. Individuals with pacemakers or internal medical devices should use caution when handling the myGaze Assistive Eye Tracker and the Mounting Bracket. Magnetic fields may affect the operation of these devices. Consult your physician and the manufacturer of your medical device to determine its susceptibility to static magnetic fields prior to handling the myGaze Assistive Eye Tracker and the Mounting Bracket. All of our magnetic products should be kept at a safe distance from individuals with these devices.

Do not handle the myGaze Assistive Eye Tracker while eating. The metal compounds in the magnets may be toxic when ingested after handling food.

## Warnings Regarding Children and Magnets

NEODYMIUM MAGNETS (RARE-EARTH) MAGNETS SHOULD BE KEPT OUT OF REACH OF CHILDREN. RARE-EARTH MAGNETS ARE NOT TOYS.

Children should not be allowed to handle or play with rare-earth magnets. Small magnets pose a choking hazard. Children and adults should not ingest magnets or place magnets in any body orifice such the ear, nose or mouth. Ingestion of magnets is very hazardous. If magnets are ingested or aspirated to the lungs, immediate medical attention is required. Swallowed magnets can stick together across intestines causing serious infections and death. Seek immediate medical attention if magnets are swallowed or inhaled. Children

under 3 should not handle magnets, in any case.

## **Transportation**

The International Air Transport Association (IATA) Dangerous Goods Regulations provide guidelines for the identification, classification, and testing of potentially hazardous materials offered for transports by air. IATA Packing Instruction 902 defines the acceptance criteria and provides packaging guidelines for magnetized material. These instructions should be consulted prior to transporting magnetic material by air. These regulations also apply to magnets built-in to products such as the myGaze Assistive Eye Tracker and the Mounting Bracket.

## 1.4 Declaration of Conformity



SMI products are for use in office environments and bear the CE mark to indicate compliance with the health and safety requirements according to European Directives. For individual product declarations please refer to [sales@smivision.de](mailto:sales@smivision.de).

All Visual Interaction eye tracking equipment has been tested and found to comply with the limits for Class B digital devices, pursuant to Part 15 of the FCC Rules and EMC directive 2004/108/EEC, and conforms to the low-voltage directive 2006/95/EEC.

## 1.5 Product Maintenance

To keep the myGaze Assistive Eye Tracker in good working order, we highly recommend that you:

- Regularly clean the shield (front glass face) of the device using the supplied microfiber cloth.
- After using the device, store it safely in the provided case.
- Do not leave the device exposed to direct sunlight, even when not in use. Store it in the case.
- Keep liquids and other contaminants away from the device.

Should the myGaze Assistive Eye Tracker become damaged, we highly recommend that you:

- Immediately unplug it from the USB port.
- Do not use the device until it has been repaired or replaced.



Do not attempt to repair the myGaze Assistive Eye Tracker by yourself. There are no user-serviceable parts in the device. Servicing, adjustment or repair should only be done by a certified distributor or by Visual Interaction GmbH.

## 1.6 Document Conventions

The following conventions are used in this document:

*Italic*

Indicates filenames and file extensions and, in some cases, product names.

**Bold**

Used for user interface buttons, selections, checkboxes, application windows and screen names.

Underlined

In the PDF and Online Help version of this manual, indicates references to a related topic in this manual or to internet addresses. In the printed version, page numbers are provided for references.



Note icon indicating additional information.



Warning icon indicating reader should pay careful attention to the information.



Reference (See Also) icon indicating a related topic. In the printed version, page numbers are shown next to topic title.

## 1.7 Glossary

### *Binocular Mode*

Both eyes are tracked to obtain values.

### *Calibration*

Ensures the gaze of the participant is accurately tracked throughout the scene by adjusting the internal eye model of the software to adapt the user's eyes.

### *Click*

In this document, it refers to a mouse click with the physical mouse. Select refers to a Gaze-Click with EyeMouse Play.

### *Calibration point*

A point or circle displayed on a screen upon which a user must fixate (or focus on) for a short duration.

### *Delay*

Sets the delay time that must elapse before a gaze dwell period starts.

### *Duration*

Sets the duration of the dwell time. This is the amount of time a user's gaze must dwell on a selection item before the item is activated.

### *Display, monitor*

Generally referring to the same device, but a Display is a term used in this document to refer to a Laptop screen, while Monitor is used in this document for a desktop PC screen.

### *Fixation*

A continuous looking at a location on a computer display.



### *Gaze*

A fixation on an element of a scene, whether on a computer display or in a real-world environment.

### *Monocular Left / Right*

Select only one eye channel. Useful when one eye cannot be tracked.

### *Operator*

Person such as teacher or caregiver who is responsible for setting up the hardware and software of an experiment using the eye tracker.

### *Operating Distance*

Distance between the subject and the cameras on the eye tracker.

### *Reference Point*

A mark on the housing of the eye tracker used to align the eye tracker with the center of the Desktop Monitor or Laptop Display.

### *Select*

In this document, it refers to a Gaze-Click with EyeMouse Play. Click refers to a mouse click with the physical mouse.

### *Tracking*

Setting which determines how the eye tracker will track the participant's eyes.

### *USB 2.0*

Generally, USB 2.0, first available in April 2000, is called High Speed USB with a limiting signaling rate of 280 Mbits/s.

# **myGaze Assistive System**

## **Chapter**




**2**

## 2. myGaze Assistive System

*myGaze Assistive System* consists of a software component, *EyeMouse Play*, and a hardware component, *myGaze® Assistive Eye Tracker*.

- *EyeMouse Play* allows users of varying abilities to access and control Windows™-based programs using only their eye gaze. With *EyeMouse Play*, tasks such as launching applications with a "double-click", selecting menu items with a "single-click" or "typing" using an on-screen keyboard can be easily performed.





- myGaze Assistive Eye Tracker tracks the user's eyes to provide input to EyeMouse Play to replicate mouse actions.




myGaze Assistive Eye Tracker mounted on a Laptop and Monitor

## 2.1 Minimum System Requirements

To achieve the best performance with myGaze Assistive Eye Tracker and EyeMouse Play, we recommend that your PC or Laptop have the following minimum system requirements:

Item	Requirements
Operating System	Microsoft™ Windows™ 7 (32 or 64 bit x86) Microsoft™ Windows™ 8 (32 or 64 bit x86) Microsoft™ Windows™ 8.1 (32 or 64 bit x86)  Do not use Microsoft XP or older.
Additional software	Microsoft DirectX Version 9.0c or later  This is installed by default in Windows. If necessary, check Microsoft Support to learn how to determine version.
PC Hardware	Intel® i5 2.6 GHz CPU or faster


Item	Requirements
	 Do not use a PC or Laptop with an AMD chip.
RAM	2 GB RAM
USB Port	USB 2.0 port
Monitor size	Screen size of between 10" to 22"



For more information see the *myGaze™ User Guide*.

## 2.2 Components

The myGaze Assistive Eye Tracker and accessories are stored in a soft carrying case to protect the device when not in use. The following table lists the components provided with the myGaze Assistive Eye Tracker.

Component		Description
myGaze Eye Tracking Device	1	 Handle the myGaze Assistive Eye Tracker with care as it is a sophisticated electronic device with sensitive cameras.
20° Angle Bracket	1	Places the myGaze Assistive Eye Tracker securely in the hinge area of a Laptop.
Mounting Bracket	1	Securely holds the myGaze Assistive Eye Tracker to the bottom side of a Monitor.
Neoprene Soft Protective Case	1	Used to safely store the myGaze Eye Tracking Device and to protect the device from scratches when stored.
Replacement magnetic tape for Angle Bracket	3	Replacement magnetic tape for using the Angle Bracket with additional Laptops.
Double-sided tape	1	Replacement tape for using the Mounting Bracket with an additional Monitor.

## 2.3 Optimal User Conditions

The myGaze Assistive Eye Tracker is an optical camera system based on infrared technology. Therefore, to ensure the myGaze Assistive Eye Tracker is operated under optimal conditions, do the following:

- The user should sit at a distance of between 50 cm and 75 cm from a Monitor or Laptop.
- Minimize any interference from direct sunlight on the myGaze Assistive Eye Tracker.
- Do not use the myGaze Assistive Eye Tracker in conditions where the user's pupils would dilate and contract frequently, such as bright lights switching on and off.
- Do not cover or block the myGaze Assistive Eye Tracker when it is powered up and is connected to the PC.
- When a user is wearing glasses, make sure the glasses are clean and free of streaks so that light does not reflect off the glasses and become visible.
- Do not use a Monitor larger than 22 inches.
- For best results, the brightness of the background color of the calibration test should be similar to the mean brightness of the stimuli shown during the experiment. This is important as to avoid large variations in the pupils of the user's eyes during the experiment, to achieve best data accuracy.

# **Getting Started**

## **Chapter**



**3**



### 3. Getting Started

Only a few easy steps are required to setup and run the myGaze Assistive System. Detailed instructions are provided in later sections of this User Guide.



Ensure your PC, Laptop, or Tablet meets the minimum system requirements. See [System Requirements](#).

#### Step 1: Obtain the Installation Package

The Installation Package contains the EyeMouse Play and required drivers.

1. Go to the **Downloads** section of the myGaze website at [forum.mygaze.com/download/](http://forum.mygaze.com/download/).
2. Enter your username and password. You can find your username and password on the invoice or delivery note from the web shop.
3. Click the download link for the EyeMouse Play installer *myGaze EyeMouse Play.msi*.
4. On the **File Download** dialog, click **Save File** to download the installer to your PC or Laptop.

#### Step 2: Install the Software

Run the EyeMouse Play installer *myGaze EyeMouse Play.msi* to start the setup wizard. Follow the instructions on the Wizard.



See [Installing myGaze EyeMouse Play](#).

## Step 2: Mount the myGaze Eye Tracker

Before using EyeMouse Play, mount the myGaze Assistive Eye Tracker to your Desktop PC or Laptop and connected to a free USB port. Carefully follow the instructions in the respective sections for a Desktop PC or Laptop.



See [Mounting the myGaze Eye Tracker](#).

## Step 4: Unhide the Gaze Positioning Icon, if required

EyeMouse Play provides a **Gaze Positioning** icon in the system tray. This may not be displayed by default. Unhide the Gaze Positioning icon, if necessary.



See [Unhiding the Gaze Positioning Icon](#)

## Step 5: Configure and Use EyeMouse Play

Start EyeMouse Play to begin using myGaze Assistive.



See [Using myGaze EyeMouse Play](#)



# **Installing EyeMouse Play**

## **Chapter**



**4**


## 4. Installing EyeMouse Play

The Installation Package for EyeMouse Play is available from the **Downloads** section of the myGaze website. This Installation Package also includes the required drivers for the myGaze Assistive Eye Tracker.

### 4.1 Downloading the Installer

To download the EyeMouse Play Installer:

1. Go to the **Downloads** section of the myGaze® website at [forum.mygaze.com/download/](http://forum.mygaze.com/download/).
2. In the **Customer Login** area, enter your log-in information provided to you by Visual Interaction GmbH.

A screenshot of a web browser displaying the 'Downloads' section of the myGaze website. The browser's address bar shows 'forum.mygaze.com/download/'. The website has a black header with the 'myGaze' logo and a navigation menu with links: HOME, PRODUCTS, SHOP, SUPPORT, LIBRARY, and COMPANY. Below the header is a large teal banner with the word 'Downloads' in white. Underneath the banner is a 'Customer Login' section. It contains a text prompt: 'Use your account data for login. You can find your username and password on the invoice or delivery note from our web shop!'. Below this is a 'Login' label and a form with two input fields: 'Username:' and 'Password:'. A 'Login' button is positioned below the password field. At the bottom right of the page, there are links for 'IMPRINT' and 'PRIVACY'.

Downloads section of the myGaze website



If you have difficulty logging in, [contact Visual Interaction](#).

3. Download the EyeMouse Play installer *myGaze EyeMouse Play.msi*.




As the installer is an executable, your browser may prevent the file from downloading. Check the browser's download settings or consult your system administrator for assistance.

## 4.2 Installing the Software

To install EyeMouse Play:


1. Run the EyeMouse Play installer *myGaze EyeMouse Play.msi*.
2. On the **Welcome** page in the Setup Wizard, click **Next** to continue.



EyeMouse Play setup wizard - Welcome


3. On the **End-User License Agreement** page, check **I agree** to accept the agreement and click **Install** to continue.
4. EyeMouse Play will now be installed.

5. When the installation has completed, click **Finish** to dismiss the Setup Wizard.



EyeMouse Play setup wizard - Installation Complete

6. A shortcut will be placed on the desktop and a **Gaze Positioning** icon will be added in the system tray.



If the **Gaze Positioning** icon did not appear in the system tray, see [Unhiding the Gaze Positioning Icon](#).

# **Mounting the myGaze Assistive Eye Tracker**

**Chapter**



**5**



## 5. Mounting the myGaze Assistive Eye Tracker

Before using EyeMouse Play, the myGaze Assistive Eye Tracker should be mounted to your PC or Laptop and connected to an available USB 2.0 port.

### Mounting on a Laptop Display

If you are using a Laptop Display, a magnetic Angle Bracket is provided to place the device in the hinge area of the Laptop.



myGaze Eye Tracker mounted in the hinge area of a Laptop



See [Placing the Device on a Desktop Monitor](#)

### Mounting on a Desktop Monitor

If you are using a PC, the Mounting Bracket is attached to the bottom edge of the Desktop Monitor.



myGaze Eye Tracker mounted on a Desktop Monitor




See [Placing the Device on a Laptop.](#)

## Mounting on a Tablet

If you are using a Tablet, you can use the magnetic Angle Bracket as well, but the device will need to be supported to sit at an angle upwards towards the user's eyes.

## 5.1 Reference Point

To ensure the device is mounted correctly, a Reference Point has been provided on the top side of the myGaze Assistive Eye Tracker. This point is used to locate the device in the center of the screen. The USB cable is the right when the device is correctly placed.



Reference Point

## 5.2 Mounting on a Desktop Monitor

The myGaze Assistive Eye Tracker is designed to be attached to the **bottom edge** (not the front face) of the Desktop Monitor frame using the Mounting Bracket. This Bracket is attached to the frame with an strong adhesive strip. Ensure the frame is clean before attaching the Mounting Bracket to ensure good adhesion.



The myGaze Assistive Eye Tracker can only be used with a Desktop Monitor with a screen size of between 10" to 22".


To mount the bracket on a Desktop Monitor:

1. Locate the horizontal center point of the Desktop Monitor (including the frame).



Mounting Bracket located at center of Monitor

2. Remove the protective cover from the adhesive strip on the bracket.



Remove adhesive strip



Replacement adhesive strips are provided in the kit for later use if you plan on using the device on other Desktop Monitors.

3. Attach the Mounting Bracket at the center point and on the bottom side of the Desktop Monitor using the adhesive strip. Press and hold for a few moments to ensure adhesion to the frame.



Mounting Bracket attached




High-strength magnets on the Mounting Bracket securely attach the myGaze Assistive Eye Tracker to the bracket. These magnets should be handled with care. See [Magnet Precautions](#).

4. Attach the myGaze Assistive Eye Tracker to the Mounting Bracket using the magnetic connectors.



Attach using magnetic connectors

5. Ensure the USB cable is to the right and the Reference Point is on top.



Reference Point on top and USB cable to the right

6. Angle the myGaze Assisive Eye Tracker upwards towards the eyes of the user.



Eye Tracker mounted on a Monitor

7. Connect the myGaze Assisive Eye Tracker to an available USB 2.0 port.
8. A completed setup of the myGaze Eye Tracker on a Desktop Monitor is illustrated below.



Completed setup

### 5.3 Mounting on a Laptop Display

An Angle Bracket is provided to use the myGaze Assistive Eye Tracker with a Laptop.

To place the myGaze Assistive Eye Tracker on a Laptop Display:

1. Locate the horizontal center point of the Laptop Display.
2. If the lower frame of the Laptop Display has sufficient height, attach a magnetic strip to the lower edge of the frame. Ensure it is centered and then connect the Angle Bracket to the magnetic strip. Otherwise, simply place the Angle Bracket in the hinge area of the Laptop.



The Angle Bracket provides a "TOP" indicator and an arrow to show the correct orientation of the bracket. This is to ensure that the Angle Bracket is located in the correct position in the hinge area of the Laptop.



Angle Bracket

3. Connect the myGaze Assistive Eye Tracker to the Angle Bracket. Ensure the USB cable is to the right and the [Reference Point](#) is on top.



Angle Bracket attached


4. Place the myGaze Assistive Eye Tracker at the horizontal center point (as located in step 1 above) and in the hinge area of the laptop.
5. Adjust the angle of the myGaze Assistive Eye Tracker by tilting the Laptop Display.





Device attached to Angle Bracket

6. Connect the myGaze Eye Tracker to the Laptop using a USB 2.0 port.
7. A completed setup of the myGaze Eye Tracker on a Laptop is illustrated below.



Completed setup



In some instances the Laptop Display may dim after connecting the myGaze Assisive Eye Tracker to the Laptop's USB port. If this happens disconnect the device from the Laptop and contact the Visual Interaction customer support line at [support@mygaze.com](mailto:support@mygaze.com).

# **EyeMouse Play Overview**


**Chapter**



**6**

## 6. EyeMouse Play Overview


EyeMouse Play consists of an **Access** button and the **myGaze Menu**. The following figure shows a typical Window™ 7 screen with EyeMouse Play running. The **Access** button is shown at the top of the screen while the **myGaze Menu** is shown in the center of the screen. The screen also dims so that the user can focus on using the **myGaze Menu**.



EyeMouse Play displaying the myGaze Menu

### 6.1 myGaze Menu

The **myGaze Menu** consists of the following functions:



myGaze Menu

- **Pause EyeMouse Play** - Stops gaze control of the mouse.
- **Gaze-Dwell Mouse Click** - Allows the user to issue a mouse-click command by fixating on an object for a pre-definable time.
- **Positioning Guide** - Provides a visual cue to the user to position themselves at the optimal position and distance from the PC Monitor or Laptop.
- **Double Click** - Allows the user to perform a double-click mouse action, such as launching an application.
- **Click and Hold** - Allows the user to click and hold an object after a dwell period.

- **Configuration Settings** - Used to configure EyeMouse Play for specific users and conditions.



The **myGaze Menu** is automatically hidden after making a selection from the menu.

## 6.2 Access Button

The **Access** button is used to show the **myGaze Menu**, which is hidden by default so as not to interfere with the use of the computer. Gazing at the **Access** button displays the **myGaze Menu**.



Access Button

Once the myGaze Assistive Eye Tracker is connected and running, and the user is sitting optimally in front of the screen, the **Access** button turns green.



Access Button - myGaze Assistive Eye Tracker Running

## 6.3 Status Colors

Three status colors are provided to help users find an optimal position to sit in front of the myGaze Assistive Eye Tracker:

- **Red** - Indicates the user is not in the correct position. Adjust the position of the user in front of the screen by using the arrows in the **Positioning Guide** window in the **myGaze Menu**.



Access button - red status

- **Orange** - Indicates the user is not in an optimal position, but can be improved. Nevertheless, the myGaze Assistive Eye Tracker is able to track the user. Adjust the position of the user by using the arrows in the **Positioning Guide** window in the **myGaze Menu**.



Access button - orange status


- **Green** - Indicates the user is in the optimal position.



Access button - green status

## 6.4 Gaze Positioning Icon

A **Gaze Positioning** icon is provided in the system tray as a status indicator. It shows the same status colors as the **Access** button.



Gaze Positioning icon in the system tray




If the **myGaze Positioning** icon does not appear in the system tray, ensure that you unhide it manually. See [Unhiding the Position Icon](#).

## 6.5 Positioning Window

EyeMouse Play provides a **Positioning** window as part of the **Calibration** window and in the **myGaze Menu** to help a user find the optimal position to




sit in front of the myGaze Assistive Eye Tracker. In the center are ovals which represent the ideal position for the user's gaze. The user simply aligns their gaze with these ovals by moving left or right. The measurement under the ovals indicates the distance (in cm) from the myGaze Assistive Eye Tracker. Arrows guide the user to move away from or closer to the myGaze Assistive Eye Tracker.



Positioning Window with user in optimum position

The color of the positioning guides, including the ovals, are based on the status colors: Red for eyes cannot be tracked, Orange for trackable eyes but not optimal and Green for trackable eyes and optimal.



Positioning Windows showing adjustment needed by the user



See [Positioning Users](#) and [Status Colors](#).

## 6.6 Overriding EyeMouse Play with PC Mouse

A care-giver or teacher can always override EyeMouse Play by simply using the PC mouse, either moving it, clicking it or right clicking it. This ensures a care-giver or teacher can always override EyeMouse Play to assist a user whenever they have difficulty performing a task. Control returns to EyeMouse Play after the PC mouse has become inactive for four (4) seconds.

## 6.7 Updating the Software

Visual Interaction GmbH relies on continuous user feedback to improve the software. Regular software updates are made available for the user through notifications in EyeMouse Play. These updates are not automatically downloaded and installed. You can decide whether you want to update the software, although it is recommended.


Notification of software updates are made available in the Title Bar of the **myGaze Menu** or on the **System** tab in the **Configuration** window.



The previous version does not need to be uninstalled.

### myGaze Menu Notification icon

When an update of EyeMouse Play is released, notification of its availability will be automatically made through a **Notification** icon in the Title Bar of the **myGaze Menu** (internet connection required). Click this icon and the update will be automatically downloaded and installed.




Notification icon for a new software release

### Systems tab notification

Alternately, if a new version is available, an update notification will appear in the **Systems** menu of the **Configuration** window:

1. Open the **Configuration** Window from the **myGaze Menu** using the **Access** button.
2. Select or click the **System** tab.
3. Information is provided about the latest version.
4. Click **Update** to install the latest version automatically.



System Update through the Configuration window

# **Using EyeMouse Play**

## **Chapter**



**7**

## 7. Using EyeMouse Play

Once the myGaze Assistive Eye Tracker has been mounted on your PC, Laptop or Tablet and EyeMouse Play has been installed, you are ready to use the myGaze Assistive System. The **myGaze Menu** provides access to the myGaze EyeMouse Play functions.




In the following sections, the term **Select** refers to a gaze, while the term **Click** refers to a mouse click with the physical mouse.

### 7.1 Unhiding the Gaze Positioning Icon

EyeMouse Play provides a **Gaze Positioning** icon located in the System Tray. This may not be displayed by default. The **Notifications** area of Windows™ is used to show the icon in the System Tray.

To show the **Gaze Positioning** icon:

1. Go to **Control Panel > Appearance and Personalization > Customize icons in the task bar**.
2. Locate **myGaze EyeMouse Play** and select **Show icon and notifications** from the dropdown menu.



Unhiding Notification icon



You can check **Always show all icons and notifications on the taskbar** so that the **Gaze Positioning** icon appears on the taskbar rather than in the System Tray pop-up window.

## 7.2 Starting EyeMouse Play


To start EyeMouse Play:

1. Ensure myGaze Assistive Eye Tracker is connected to a USB 2.0 port.
2. Start EyeMouse Play.



The screen dimensions of the Laptop, PC or Tablet will be automatically detected.


3. The **Access** button and the **myGaze Menu** will appear with the desktop dimmed to allow the user to focus on the **myGaze Menu**.



myGaze Assistive Eye Tracker starting

1. Once the myGaze Assistive Eye Tracker has started, use the **Positioning** area in the **myGaze Menu** to position the user.





myGaze Menu and Access Button displayed



See [Positioning Users](#)

4. Perform a calibration (recommended). Otherwise a 0-point calibration will be done automatically.



You can configure myGaze EyeMouse Play to always calibrate on

startup. See [Setting Calibration Properties](#).



See [Performing a Calibration](#).

5. The user can now select an item using only their eye gaze.
6. After making a selection, the **myGaze Menu** will be automatically hidden.



A user profile can be set so that on subsequent launches of EyeMouse Play, the last saved calibration will be used instead of performing another calibration. See [Managing Profiles](#).


## 7.3 Positioning Users

A user needs to be positioned in front of the myGaze Assistive Eye Tracker so the device can reliably track their eyes.

Here are a few guidelines for positioning users:


- Ensure the myGaze Assistive Eye Tracker is oriented upwards towards the user's eyes. It should be directed towards the user's eyes as the user sits comfortably in front of the Desktop Monitor or Laptop Display.
- Users should sit approximately 50 cm to 75 cm from the screen (optimal user position is between 60 cm and 70 cm away).
- Use the **Positioning Guide** icon in the system tray. When its color is green, the user's gaze can be tracked.

The following shows a user sitting centered and approximately 51 cm away from the myGaze Assistive Eye Tracker.



Positioning Window - Green (optimal)

The following shows a user sitting off center to the left and too close to the screen. Arrows are provided to guide the user to move to the right and away from the screen.



Move to the right and away from the screen

## 7.4 Opening the myGaze Menu

The **myGaze Menu** can be opened in several ways.

- **Access button** - A user gazes at the **Access** button.
- **Gaze point** - A user gazes at the center area of the myGaze Assistive Eye Tracker.
- **Gaze Positioning Icon** - A caregiver or teacher can open the **myGaze Menu** using the popup menu from the **Gaze Positioning** icon in the System Tray.

### Access Button

By gazing at the **Access** button, the **myGaze Menu** opens. The screen is dimmed to allow the user to focus on the **myGaze Menu**.




Gaze at Access button to open  
myGaze Menu



**Access** button properties, including the gaze duration, are set in the **Settings** tab of the **Configuration** window. See [Setting Access Properties](#).

### Gaze Point


By gazing at the center point of the myGaze Assistive Eye Tracker, the **myGaze Menu** opens. The screen is dimmed to allow the user to focus on the **myGaze Menu**.



Gaze point on the Assitive Eye Tracker

## System Tray

Right click on the **Gaze Positioning** icon in the System Tray. In the popup menu, select **Open** to display the **myGaze Menu**.



Accessing from Status Bar



If the **Positioning Guide** icon is not shown in the system tray, see [Unhiding the Gaze Positioning Icon](#).

## 7.5 Selecting an Item from myGaze Menu

To select an item from the **myGaze Menu**:

1. Gaze at the **Access** button to open the **myGaze Menu**.



Access Button

2. Gaze at an item in the **myGaze Menu**.
3. When the item is selected, a checkbox will appear in the selected icon and the **myGaze Menu** will be automatically hidden.


## 7.6 Activating Cursor Mode

When EyeMouse Play is running, the **Cursor** mode is used. In this mode, the user controls the movements of the cursor using only their eye gaze.

To switch between another mode and **Cursor** mode:

1. Gaze at the **Access** button to open the **myGaze Menu**.

2. Deselect the currently used mode by gazing at the respective button.
3. At the end of the Dwell-Gaze duration, the checkmark in the respective button disappears indicating **Cursor** mode is now active.



myGaze EyeMouse Play in Cursor mode

5. The **myGaze Menu** will be automatically hidden.


## 7.7 Pausing EyeMouse Play

Control of the mouse pointer can be returned from EyeMouse Play to the physical mouse by pausing EyeMouse Play. When activated, the **Paused** menu item is unchecked.

To pause EyeMouse Play:

1. Gaze at the **Access** button to open the **myGaze Menu**.
2. Gaze at the **Pause EyeMouse Play** button.
3. At the end of the Dwell-Gaze duration, the checkmark in the **Pause EyeMouse Play** disappears indicating **EyeMouse Play** is paused.






4. The **myGaze Menu** will be automatically hidden.

## 7.8 Activating Dwell-Gaze Mouse Click Mode

A single-click mouse command is enabled by the **Dwell-Gaze Mouse Click** button on the **myGaze Menu**. **Dwell-Gaze Mouse Click** mode allows a user to perform tasks that require this, such as clicking a menu item. The fixation time for gazing at an item on the screen is called the **Dwell-Duration** and is a pre-defined time set in the **Configuration** window.

To activate **Dwell-Gaze Mouse Click**:

1. Gaze at the **Access** button to open the **myGaze Menu**.
2. Gaze at the **Dwell-Gaze Mouse Click** button.
3. At the end of the Dwell-Gaze duration, a checkmark in the **Dwell-Gaze Mouse Click** button appears indicating **Mouse Click Mode** is active.



Dwell-Gaze Mouse Click activated

3. The **myGaze Menu** will be automatically hidden.




To set the dwell-gaze duration, see [Setting Mouse Gaze Functions](#)

## 7.9 Activating Double-Click Mode

A double-click mouse command is enabled by the **Double-Click** button on the **myGaze Menu**. **Double-Click** mode allows the user to perform tasks that requires this, such as starting a program from a Desktop icon.

To activate **Double-Click** mode:

1. Gaze at the **Access** button to open the **myGaze Menu**.
2. Gaze at the **Double-Click** button.
3. At the end of the Dwell-Gaze duration, a checkmark in the **Double-Click** button appears indicating **Double-Click Mode** is active.




3. The **myGaze Menu** will be automatically hidden.

## 7.10 Activating Click and Hold Mode

A Click and Hold mouse command is enabled by the **Click and Hold** button on the **myGaze Menu**. **Click and Hold** mode allows the user to click and hold on an object to move it around the screen or application. The object moves with the user's gaze. For example, in a paint program, a user could paint by moving a "paint brush" around with their eyes to draw a picture.

To activate **Click and Hold** mode:

1. Gaze at the **Access** button to open the **myGaze Menu**.
2. Gaze at the **Click and Hold** button.
3. At the end of the Dwell-Gaze duration, a checkmark in the **Click and Hold** button appears indicating **Click and Hold** is active.




Activating Click and Hold mode

3. The **myGaze Menu** will be automatically hidden.

## 7.11 Shutting Down EyeMouse Play

To shut down EyeMouse Play:

1. Using the PC Mouse, right-click on the **Gaze Positioning** icon in the system tray.



Select Exit to shut down myGaze  
EyeMouse Play

2. Click **Exit**.

# **Configuring EyeMouse Play**

## **Chapter**




**8**



## 8. Configuring EyeMouse Play

Teachers or caregivers can use the **Configuration** window to set up a profile and configure EyeMouse Play for specific users and conditions.



Configuration window

The **myGaze Menu** consists of the following functions:

- **Manage profiles** - Creates a user profile so that on subsequent launches of EyeMouse Play, the last saved calibration and other settings will be used. See [Managing Profiles](#).
- **Perform calibrations** - Adapts the internal software of the myGaze Assistive Eye Tracker to the unique characteristics of a user's eyes. See [Performing a Calibration](#).
- **Setting EyeMouse Play gaze functions** - Configures the mouse gaze to


suit the needs and abilities of the user. See [Configuring Mouse Settings](#).

- **Setting keyboard shortcuts** - Provides settings for keyboard shortcuts for EyeMouse Play functions. See [Setting Keyboard Shortcuts](#).
- **System Settings** - Provides a number of system related settings for EyeMouse Play. See [System Settings](#).

## 8.1 Viewing Configuration Window

To view the **Configuration** window:

1. Click or gaze at the **Configuration Settings** button.



2. The **Configuration** window appears.

## 8.2 Managing Profiles


A Default user profile is provided that includes a number of settings which are suitable for most users.

However, you can create a user profile that can be used on subsequent launches of EyeMouse Play. This profile will be use each time EyeMouse Play is run until another profile is selected. This can be useful if you want to use the last saved calibration instead of performing another calibration and if you want to use the same settings.

### 8.2.1 Creating a Profile


To create a user profile:

1. In the **myGaze Menu**, click the **Configuration Settings** button to open the **Configuration** window.
2. Click the **Default** button in the **Configuration** window to open the **Select profile** dialog.



Select Profile dialog

3. Click **Add** to open the **Add new profile** window




Add a new profile

4. Enter a new profile name and click **OK**.
5. You can now modify the settings for this new profile.

### 8.2.2 Selecting a Profile

To switch between profiles:

1. Click the currently set profile to open the **Select a Profile** dialog.



Select from the available profiles

2. Select an available profile to load the settings of this profile.


### 8.2.3 Deleting a Profile



A currently loaded profile cannot be deleted. You need to change to another profile before deleting a profile.

To delete a profile:

1. Change to another profile. See [Selecting a Profile](#).
2. Click the **Delete** button to open the **Select profile** dialog.



Delete a selected profile

3. All available profiles that can be deleted will be shown.
4. Click the profile you want to delete
5. Confirm your deletion.
6. All the settings and profile name will be deleted.

## 8.3 Performing a Calibration


Calibration adapts the internal software of the myGaze Assistive Eye Tracker to the unique characteristics of a user's eyes. This is done to achieve the best possible data accuracy. A successful calibration ensures that the myGaze Assistive Eye Tracker accurately tracks where the user is looking on the screen.

To perform a calibration, a user fixates on a series of targets that are sequentially displayed on the screen.

### 8.3.1 Calibration Methods

myGaze EyeMouse Play provides five different calibration methods with increasing accuracy as more fixation points are provided:


- **0-Point Calibration** (Calibration-free mode) - This mode uses pre-set data. This mode can be used for most users.
- **1-Point Calibration** - This mode uses a single fixation point at which the user focuses on to calibrate the eye gaze.



1-point calibration


- **2-Point Calibration** - This mode uses two fixation points that the user focuses on in succession to calibrate the eye gaze.





2-point calibration

- **5-Point Calibration** - This mode uses a set of five fixation points that the user focuses on in succession to calibrate the eye gaze.



5-point calibration


- **9-Point Calibration** - This mode uses a set of nine fixation points and

provides optimal accuracy.

### 8.3.2 Cascading Points

To help the user to locate and fixate on a calibration point, the point cascades from a large circle to a small circle. When the user has correctly fixated on the point, the calibration will stop when the point is at its smallest. The next point will then be displayed.

The following shows a typical cascading point.




Cascading points

### 8.3.3 Setting Calibration Properties

The **Calibration** tab provides settings for a calibration that best suits the abilities of a user. This includes setting the calibration method and the properties of the calibration fixation points at which a specific user must gaze for the required gaze time.

To set **Calibration** properties:

1. In the **myGaze Menu**, click the **Configuration Settings** button to open the **Configuration** window.
2. Select the **Calibration** tab.




Calibration tab

3. Set the properties as required.
4. Click or select **OK** to dismiss the **Configuration** window.


**Calibration** tab settings include:

- **Calibration Method** - Select from **0 Point**, **1 Point**, **2 Point**, **5 Point** and **9 Point**.




Selecting a calibration method

- **Calibration Point** - Select from **Circle**, **Crosshair**, **Star**, or an **Image**. By selecting **Image**, a Windows **Browse** dialog appears to load a suitable image for the point.




Select a calibration point

- **Color** - Select from seven predefined colors.



Select a calibration point color

- **Sound** - Check to play a ping sound when a fixation point has been accepted. Default is **Play a sound**.
- **Tracking** - Select from **Binocular** (both eyes are tracked), **Monocular Left** (only left eye is tracked) or **Monocular Right** (only right eye is tracked).



Select tracking mode

- **Accept Points Automatically** - Sets acceptance method of calibration points. If set to **Automatic** (selected), the point is accepted when the software identifies that the user's gaze has fixated on a point at the moment when it is at its smallest size. If set to **Manual** (deselected), the user presses the **Spacebar** to accept a point.
- **Always calibrate at startup** - Causes the **Calibration** window to appear when starting EyeMouse Play.

### 8.3.4 Running a Calibration

Calibration is performed in the **Calibration** window in full screen mode. This window is accessed through the **Calibrate** button in the **Configuration** window. After selecting a Calibration method, a series of calibration points are displayed in succession, depending on the calibration method selected. The gaze of a calibration point is accepted either automatically or manually. If set

to manually, the user presses the Spacebar to accept a calibration point.


Acceptance normally occurs when the user fixates on a calibration point long enough for the software to determine the gaze location of a user. If a user has difficulty fixating on a point, using manual acceptance is preferred so the user can take sufficient time as required to fixate on a given point.



Running a calibration can be set to **Always calibrate on startup**. See [Set Calibration Properties](#).

To run a calibration:

1. From the **Configuration** window, select the **Calibration** tab.




Calibration tab

2. Click **Calibrate**.



3. EyeMouse Play goes into full screen mode displaying the **Calibration** window.




4. Position the user so that the white ovals are within green outline ovals. This indicates the user is sitting in the correct position in front of the myGaze Assistive Eye Tracker.



Visual cues, including arrows, colors and distance measurements are provided to help the user sit in the optimal position. see [Positioning Users Correctly](#).

5. Select a Calibration method: 1 point, 2 point, 5 point or 9 point.
6. Click **Calibrate** to begin the calibration.
7. After several seconds, fixation points will be displayed in succession on the screen, depending on which calibration mode was selected. The user must focus on each point as it is displayed. The following shows five fixation points for a 2-point calibration.




5 Point Calibration

8. If **Accept Points Manually** was set in the **Configuration** window (that is, **Accept points automatically** was deselected), click the **Spacebar** to accept a point when the circle is at the smallest size. If **Accept points automatically** was set (the default setting), acceptance is done when the software determines that the user has correctly fixated on the point.



See [Set Calibration Properties](#).

9. The results of the calibration will be shown.



Calibration results

10. If the results are unusual or inadequate, click **Improve** to run the calibration test again on the failed calibration points, while ensuring the user keeps a focus on each point as it appears. See [Improving a Calibration](#).


11. After calibration has been completed, click **Finish**.

### 8.3.5 Improving a Calibration

During the calibration process, if the system encounters 3 failed attempts to calibrate on a point, this point will be disregarded and move on to the next point.

At the end of the process the user can see how many points have been successfully calibrated and choose to improve their calibration by recalibrating the omitted points.

In the following example, 2 Point Calibration was used, as indicated by the two stars in the upper left corner. However, one point failed to calibrate properly. By clicking **Improve**, the user can recalibrate on the failed point.




## 8.4 Setting Mouse Properties

The **Mouse** tab provides settings to configure the mouse gaze to suit the needs and abilities of the user.

To set **Mouse Play** properties:

1. In the **myGaze Menu**, click the **Configuration Settings** button to open the **Configuration** window.
2. Select the **Mouse** tab.



Mouse tab settings

3. Set the properties as required.
4. Click or select **OK** to dismiss the **Configuration** window.




A live feedback of the changes is shown to verify the settings.



Live Feedback window

**Mouse** tab settings include:

- **Smooth** - Assists the user in gazing around a screen. If smoothing is set to a lower value, then the gaze cursor will have a more abrupt movement across the screen. A higher smoothing value will make the gaze cursor have a more stable, less jumpy movement.
- **Activation method** - Sets which method causes an action to occur. **Dwell** causes a mouse action to occur after the set dwell time. **Blink** causes a mouse action to occur after the user blinks. **Switch** causes a mouse action to occur after a user hits a switch, which is a set keyboard shortcut (when **Switch** is selected, an option appears to set a keyboard shortcut).



Select Activation Method

- **Size** -Sets the size of the dwell circle. Settings must be between 10 px and 160 px.
- **Delay** - Sets the delay time that must elapse before a gaze dwell period starts. For example, if the delay time is set to 2 seconds, then the dwell duration will only begin after the user has gazed at a button or menu item

for 2 seconds.


- **Duration** - Sets the duration of the dwell time. It can be increased or decreased from a quick pulse of 0.4 seconds to a longer duration of 2.0 seconds. This is the amount of time a user's gaze must dwell on a selection item before the item is activated. For example, if **Duration** is set to 1 second, then a user must gazes at a button or menu item for 1 second before the application responds to the user request.

## 8.5 Setting Access Properties

The **Access** tab provides settings for how the access the **myGaze Menu**. This can be through either the **Access** button or through the Gaze point on the myGaze Assistive Eye Tracker. The **Access** button can also be disabled and hidden.

To set **Access** properties:

1. In the **myGaze Menu**, click the **Configuration Settings** button to open the **Configuration** window.
2. Select the **Access** tab.




Access tab

3. Set the properties as required.
4. Click or select **OK** to dismiss the **Configuration** window.

**Access** tab settings include:

- **Gaze at camera** - Set to use the gaze point of the myGaze Assistive Eye Tracker to show or hide the **myGaze Menu**. Deselect to use the **Access** button.
- **Access Button** - Select from the dropdown list the most convenient location on the screen for the **Access** button. Selections include **Off, Top, Bottom, Left, Right**. If you do not want to use the **Access** button but only the gaze point of the myGaze Assistive Eye Tracker, select **Off** to hide the **Access** button.





Access button position settings

- **Size** - Set the size of the **Access** button. Use the slider to change the size.

## 8.6 Setting Keyboard Shortcuts

The **Shortcuts** tab provides settings for keyboard shortcuts for EyeMouse Play functions. With a keyboard shortcut a caregiver or teacher can manage EyeMouse Play functions without interrupting the program or activities of the user.


For instance, dwell time can be reduced or increased by using a key at the same time as the user continues with a task uninterrupted.

To set keyboard shortcuts:

1. In the **myGaze Menu**, click the **Configuration Settings** button to open

the **Configuration** window.

2. Select the **Shortcuts** tab.



Keyboard Shortcuts tab

3. Click the dropdown menu to the right of an item to select an option.



**None** indicates that no keyboard shortcut has been assigned. If a key has been assigned, this key will appear in the dropdown menu.

4. Press any key to set the shortcut.
5. Click or select **OK** to dismiss the **Configuration** window.



Be careful which key is chosen as a keyboard shortcut. Some keys


may be reserved for some applications and when pressed during the use of such an application can cause unexpected behaviors. Check for reserved keys for specific applications before setting.

## 8.7 System Settings

The **System** tab provides a number of system related settings for EyeMouse Play.

To set **System** properties:

1. In the **myGaze Menu**, click the **Configuration Settings** button to open the **Configuration** window.
2. Select the **System** tab.



System tab

3. Set the properties as required.
4. Click or select **OK** to dismiss the **Configuration** window.

**System** tab settings include:

- **Launch at Windows startup** - Select to start EyeMouse Play after Windows has started and the myGaze Assistive Eye Tracker has been powered up and running.
- **Display** - If you are using two monitors, you can select either Primary or Secondary monitor and EyeMouse Play will automatically set the geometry values.
- **Powersave** - When selected, the Powersave mode allows the system to 'sleep' when the user is not in front to the device. When the user's face is recognized, the system becomes active again after a few moments. This mode is best suited for tablets and when the device is running on battery power.

# **EyeMouse Play and Third-Party Software**

**Chapter**



**9**

## 9. EyeMouse Play and Third-Party Software

myGaze Assistive works well with a variety of third-party assistive programs for curriculum study, basic learning, AAC (Augmentative and Alternative Communication), leisure and rehabilitation activities, and more. This section describes Grid 2 AAC software, a popular AAC program by Sensory Software.



More information about suitable third party programs can be found on [www.myGaze.com](http://www.myGaze.com) as well as from our global reseller network.

### 9.1 Grid 2 AAC Overlay

myGaze Assistive works out-of-the box with *The Grid 2* software, a popular AAC (Augmentative and Alternative Communication) program by Sensory Software. With this program, a user can rely on EyeMouse Play to access Grid 2 functions through their eye gaze.

#### Using Grid 2 with EyeMouse Play

EyeMouse Play has only two modes within the Grid 2:

- **Paused mode** - cannot issue any click commands to Grid 2.



Paused Mode

- **Active mode** - any one of [Dwell-Gaze Mouse Click](#), [Double-Click](#), [Click and Hold](#) as well as the [default Cursor mode](#) which issues a Dwell-Click to Grid 2.

To use Grid 2 AAC Overlays with EyeMouse Play:

1. Switch to Active mode by selecting any of the four Active commands in the **myGaze Menu**. In each of those modes EyeMouse Play will issue a Dwell-Click command within the Grid 2.
2. Select a button on the on-screen keyboard by looking at the button.
3. When the **Dwell-Gaze** period is complete, a click will occur.




Dwell-Click in Grid

## Independently Accessing myGaze Menu in The Grid 2

Accessing **myGaze Menu** in Grid 2 can be done through the two standard options:

- **Gaze Point** - Gazing at Gaze point of the myGaze Assistive Eye Tracker opens the myGaze menu is Grid 2.



Gaze point

- **Access** button - Gazing at the **Access** button opens the **myGaze Menu** in Grid 2.




Access button



If you choose the **Access** button option, set the size and location of the **Access** button so it will not interfere with the onscreen grid. See [Setting Access Button Properties](#).





Access button shown overlaying the bottom part of The Grid 2 application


## Customizing the Landing Grid

In addition to the standard access options, a user can also easily customize the landing grid for EyeMouse Play. This is done by adding grid buttons for all or some of the EyeMouse Play commands. These commands include:

- Calibration
- Positioning guide
- Gaze menu
- Configuration



A standard myGaze grid created by Visual Interaction is available for installation in the EyeMouse Play folder after EyeMouse Play has been installed.



myGaze Grid

# **Appendix**

## **Chapter**



**10**

## 10. Appendix

This section provides troubleshooting information as well as license agreements.

### 10.1 Troubleshooting

#### *Time required to start up and configure EyeMouse Play*

Please be aware that the EyeMouse Play needs time to start up and configure the device. The Eye Tracking system is not ready until it enables the IR illumination and starts to track the participant.

#### *Restarting EyeMouse Play*

It can be helpful in certain cases to simply quit the EyeMouse Play application and restart it again.

#### *Dual monitor issues*

The EyeMouse Play application is able to use a dual monitor setup. Be sure that the calibration of the Eye Tracking device will be performed on the same monitor which will be used by the user.

#### *myGaze Eye Tracking Software issues.*

For any questions regarding the myGaze Eye Tracking software please read the manuals delivered with the corresponding software.

## 10.2 License Agreement and Warranty

### License Agreement and Warranty for Visual Interaction Software Products

DO NOT CLICK "I AGREE" UNLESS (1) YOU ARE AUTHORIZED TO AGREE TO THE TERMS AND CONDITIONS OF THIS LICENSE AGREEMENT ON BEHALF OF YOURSELF AND/OR YOUR COMPANY; AND (2) YOU INTEND TO ENTER INTO THIS LEGALLY BINDING AGREEMENT ON BEHALF OF YOURSELF AND/OR YOUR COMPANY.

**IMPORTANT – PLEASE READ CAREFULLY:** This license agreement ("Agreement") is an agreement between you (either an individual, an institution or a company, "Licensee") and Visual Interaction GmbH ("VI"). The "Licensed Materials" provided to Licensee subject to this Agreement include the software product purchased by you ("SOFTWARE") as well as any "on-line" or electronic documentation associated with the SOFTWARE, or any portion thereof ("Documentation"), as well as any updates or upgrades to the SOFTWARE and Documentation, if any, or any portion thereof, provided to Licensee at VI's sole discretion.

By installing, downloading, copying or otherwise using the Licensed Materials, you agree to abide by the following provisions. This Agreement is displayed for you to read prior to using the Licensed Materials.

If you do not agree with these provisions, do not download, install or use the Licensed Materials.

1. License. Subject to the terms of this Agreement, VI hereby grants and Licensee accepts a non-exclusive license without the right to sublicense ("License") to use the Licensed Materials only in accordance with the Documentation provided by VI. License is subject to Licensee providing to or allowing VI to collect from time to time certain data that is collected by the Licensed Materials ("Performance Data") for purposes of quality control and improvement. Installation of the SOFTWARE is Licensee's sole responsibility.
2. Rights in Licensed Materials. Unless otherwise specified in the product description, this License does not grant any rights to obtaining future

upgrades, updates or supplements of the Software. If upgrades, updates or supplements of the Software are obtained, however, the use of such upgrades or updates is governed by this License Agreement and the amendments that may accompany them and may be subject to additional payments and conditions.

VI retains all right title and interest in and to the SOFTWARE and to Performance Data, and no other license is granted to Licensee by implication, estoppel or otherwise.

In particular and without prejudice to mandatory rights under statutory law, Licensee agrees not to:

- (i) copy or use the SOFTWARE and Performance Data in any manner except as expressly permitted in this Agreement, including – but not limited to – duplicating a single user License;
- (ii) attempt to circumvent, disable or defeat the limitations on Licensee's use of SOFTWARE,
- (iii) rent, lease, distribute, sublicense or loan the SOFTWARE in whole or in part to any third party or make the SOFTWARE publicly available;
- (iv) use the SOFTWARE for providing time-sharing services, any software-as-a-service offering ("SaaS"), service bureau services or as part of an application services provider or as a service offering; or
- (v) alter or remove any proprietary notices in the SOFTWARE.

3. Confidentiality. Licensed Materials and Performance Data are proprietary to VI and constitute VI trade secrets. Licensee shall maintain Licensed Materials and Performance Data in confidence and prevent their disclosure using at least the same degree of care it uses for its own trade secrets, but in no event less than a reasonable degree of care. Licensee shall not disclose Licensed Materials and Performance Data or any part thereof to anyone for any purpose, other than to its employees and sub-contractors for the purpose of exercising the rights expressly granted under this Agreement, provided they have in writing agreed to confidentiality obligations at least equivalent to the obligations stated herein. Confidentiality obligations do not apply if and insofar

- (i) the information became available to Licensee on a non-confidential basis from a source other than VI or its licensors or affiliates, provided that the source is entitled to disclose the information or
  - (ii) the information was or became generally available to the public other than as a result of a disclosure by a licensee (including the Licensee) in breach of his contractual obligations.
4. Limited Warranty and Liability. The provisions of VI's current version of the General Terms and Conditions for the Supply of Goods and Services apply with regard to VI's limitation of warranty and liability. These limitations of liability shall also apply in case of Licensee's claims for damages against VI's employees or agents.
- Pursuant and subject to the General Terms and Conditions for the Supply of Goods and Services VI only warrants that the SOFTWARE is free from defects and/or third party rights (as specified in those terms).
5. Licensee Indemnity. Licensee will defend and indemnify VI, and hold it harmless from all costs, including attorney's fees, arising from any claim that may be made against VI by any third party as a result of Licensee's use of Licensed Materials.
6. Export Restriction. Licensee will not remove or export from Germany or from the country Licensed Materials which were originally transmitted and/or shipped to by VI or re-export from anywhere any part of the Licensed Materials or any direct product of the SOFTWARE except in compliance with all applicable export laws and regulations, including without limitation, those of the U.S. Department of Commerce.
7. Non-Waiver; Severability; Non-Assignment. The delay or failure of either party to exercise any right provided in this Agreement shall not be deemed a waiver.

If any provision of this Agreement is held invalid, all others shall remain in force. Licensee may not, in whole or in part, assign or otherwise transfer this Agreement or any of its rights or obligations hereunder.

8. Entire Agreement; Written Form Requirement. There are no separate oral agreements; any supplementary agreements or modifications to this License Agreement must be made in writing (whereas for the General Terms and Conditions for the Supply of Goods and Services and the purchase process electronic form is sufficient). This also applies to any waiver of this requirement of written form.
9. Notices. All notices under the Agreement must be in writing and shall be delivered by hand or by overnight courier to the addresses of the parties set forth above.
10. Applicable Law and Jurisdiction. German law with the exception of its conflict of laws rules shall apply exclusively. The application of the United Nations Convention on Contracts for the International Sale of Goods (CISG) is expressly excluded. The courts of Potsdam, Germany, shall have exclusive jurisdiction for any action brought under or in connection with this Agreement.

© 2014 Visual Interaction GmbH  
Potsdam, Germany

## 10.3 Third Party Licenses

This software is based on pugixml library (<http://pugixml.org>).  
Copyright (C) 2006-2014 Arseny Kapoulkine.

This software is based on zlib/libpng library (<http://www.libpng.org/> / <http://www.zlib.net/>).  
Copyright (C) 1995-2013 Jean-loup Gailly and Mark Adler



# About Visual Interaction

## Chapter



11

## 11. About Visual Interaction

Visual Interaction (VI) makes gaze-based interaction, multimodal user interfaces and focused analysis of visual interaction commonplace and affordable. Based on leading technology from SMI, a leader in eye tracking for 20 years, Visual Interaction brings to market myGaze, an easy to use and cost efficient stationary gaze tracking solution, specifically customized for the need of interactive solutions. For more information, see [www.myGaze.com](http://www.myGaze.com).

# Index

## - 0 -

0-Point Calibration 76

## - 1 -

1-Point Calibration 76

## - 2 -

2-Point Calibration 76

## - 5 -

5-Point Calibration 76

## - 9 -

9-Point Calibration 76

## - A -

Accept gaze points

Automatically 78

Manually 78

Access button

Location of 91

Opening myGaze Menu 56

Size of 91

Access settings 91

Accessing myGaze Menu 56

Activation Method

Blink 88

Dwell 88

Switch 88

Always Calibrate at startup 78

Arrows

Use of in Positioning Guide 54

Assisting users to position correctly 54

Assistive Eye Tracker

Position of user in front of 54

Automatically accepting gaze points 78

## - B -

Binocular Mode 11

Binocular tracking 78

Blink

Activation method 88

## - C -

Calibrate

Use last saved 78

Calibrating

Always on startup 78

Calibration 78

About 75

Defined 11

Improving 87

Methods 76

Recalibrating a point 87

Running 83

Setting up 78

Sound used in 78

Calibration method

Selecting 78

Calibration point 78

Defined 11

Calibration Points

Cascading 78

Calibration-free mode 76

Cascading calibration points 78

Click and Hold 65

Color of fixation point 78

Components

## Components

myGaze Assistive System 18

## Configuration settings

Access 91

keyboard shortcuts 93

Mouse 88

System 95

## Configuration window

Viewing 70

## Configuring

EyeMouse Play 69

## Contacting Visual Interaction 3

correct position in front of Eye Tracker 54

CPU 16

## Cursor

Switching to 58

## - D -

Declaration of Conformity 8

Delay 11, 88

Dilating pupils 18

Direct sunlight, avoid 18

DirectX 16

## Display

Selecting 95

## Distance

Between user and Eye Tracker 18

Distance from Eye Tracker 54

Document conventions 10

Document version 2

Double-Click button 64

## Downloading

EyeMouse Play installer 25

Dual Monitors 104

Duration 11

## Dwell

Activation method 88

Dwell Duration/Time 88

Dwell-Gaze Mouse Click button 62

Dwell-time

Set in Configuration Menu 62

## - E -

## EyeMouse Play

Configuring 69

Grid 2 AAC Overlay 98

Installing 26

Overriding with PC mouse 46

Overview 40

Restarting 104

Shutting down 67

Starting 51

Third-Party Software 98

Updating 46

Using 50

## - F -

## Fixation

Defined 11

## - G -

Gaze Dwell Size 88

## Gaze

Defined 11

Gaze at camera setting 91

Gaze control

Stopping 60

Gaze Point

Opening myGaze Menu 56

Gaze Points

Accepting 78

- 
- Gaze Positioning icon 44
    - Exit EyeMouse Play 67
    - Unhiding 50
  - Glasses
    - User wearing 18
  - Glossary 11
  - Grid 2 AAC Overlay 98
  - H -**
  - How to set up myGaze Assistive System 21
  - I -**
  - Icons
    - Explained 10
  - Installer
    - EyeMouse Play obtaining 25
  - Installing EyeMouse Play 25
  - Intel® i5 16
  - K -**
  - keyboard shortcut
    - Calibration 93
    - Double click 93
    - Dwell duration - 93
    - Dwell duration + 93
    - Left click 93
    - Pause 93
    - Positioning guide 93
    - setting 93
  - L -**
  - Laptop
    - myGaze Assistive Eye Tracker 35
  - Launch at Windows startup 95
  - Launching myGaze EyeMouse Play at Windows startup 95
  - License Agreement and Warranty 105
  - Licenses, pugixml, zlib, libpng 108
  - Live Feedback window 88
  - Location of Access button 91
  - M -**
  - Magenets
    - Warnings about 5
  - Manually accepting gaze points 78
  - Microsoft® DirectX 16
  - Monitor
    - Maximum size of 18
    - Minimum size of 18
    - myGaze Assistive Eye Tracker 31
  - Monitors
    - Primary or secondary selecting 95
  - Monocular Left / Right
    - Defined 11
  - Monocular left tracking 78
  - Monocular right tracking 78
  - Mounting 29
  - Mounting Bracker 31
  - Mouse button 88
  - Mouse settings
    - Configuring 88
  - Mouse tab 88
  - Mouse-click command 62
  - myGaze 2
  - myGaze Assistive
    - Getting started 21
  - myGaze Assistive Eye Tracker 29
    - and opening myGaze Menu 56
  - Liability Information 4
  - Magnets 5
  - Maintaining 9

myGaze Assistive Eye Tracker 29

Mounting 29

Mounting on Laptop 35

Mounting on monitor 31

Reference point 30

Use conditions 18

myGaze Assistive System

Components 18

Overview 15

myGaze EyeMouse Play

Launching at startup 95

myGaze Menu 42

Access button 42

Click and Hold 65

Configuration Settings button 70

Double-Click 64

Dwell-Gaze Mouse Click 62

Hiding 42

Opening 56

Overview 40

Pausing EyeMouse Play 60

Selecting item in 58

Showing 42

Status colors 43

Switching to Cursor 58

myGaze website 2

**- N -**

Neodymium magnets

Safety information about 5

Note icon 10

**- O -**

Opening myGaze Menu 56

Using Access button 56

Using Gaze Point 56

opening myGaze Menu, using Positioning Guide 56

Operating Distance 11

Operating System 16

Operator 11

Optimal user conditions 18

Overriding EyeMouse Play with PC mouse 46

**- P -**

Pausing EyeMouse Play 60

PC Monitor

Maximum size of 18

PC mouse

Overriding EyeMouse Play 46

Returning control to mouse 60

Positioning Guide, to open myGaze Menu 56

Positioning Window

About 44

PowerSave mode 95

Profile

Creating 71

Profiles

Deleting 74

Managing 71

Selecting 73

**- R -**

RAM 16

Recalibrating 87

Reference point 30

Defined 11

Release date 2

Requirements

System 16

---

Restarting EyeMouse Play 104

## - S -

See also icon 10

Setting keyboard shortcuts 93

Smoothing 88

Software version 2

Sound

Used in calibration 78

Status colors

Green 43

Orange 43

Red 43

Switch

Activation method 88

System

Requirements 16

System tray

Gaze Position icon showing in 50

## - T -

Technical Support 3

Third Party Licenses 108

Third party software 98

Toggling between Cursor and Gaze-Dwell  
Mouse Click 58

Tracking

Defined 11

tracking, binocular 78

tracking, in calibration 78

tracking, monocular left 78

tracking, monocular right 78

Troubleshooting 104

Typographical conventions 10

## - U -

Unhiding Gaze Positioning icon 50

USB 16

USB 2.0

Defined 11

Requirements 16

User conditions

Optimal 18

Using dual monitors 104

## - V -

Visual Interaction

About 110

Contacting 3

## - W -

Warning icon 10

Wearing glasses 18

Windows® 7 16

Windows® 8 16

Windows® XP 16

Wizard

EyeMouse Play setup 25

